


## QuickONE PAD

mobile 2D barcode scanner  
and data collector



QuickONE PAD is the first mobile data collector for 2D barcodes with an integrated optical unit. It is possible to read single tubes with data matrix barcodes on the bottom and collect their data for later use. The enhanced export module allows directly generating reports with the captured codes and their individual timestamps for documentation.

### QuickONE PAD Components

#### Data Collection

- mobile 2D barcode data collector
- optional sample to patient assignment

#### Standard Export

- export files in csv, pdf and html format

#### User Interface

- Intuitive operation with touch optimized interface
- sample highlighting (pos / neg)

#### mobile 2D barcode reader

The attached 2D Barcode Reader allows direct scanning of barcodes to the tablet.

All key functions will be displayed around the reader module to provide a user interface optimized for touch input.



#### data collector


Captured barcodes are listed with their individual timestamps. With the enhanced data collection module, it is also possible to collect additional patient information (like patient ID, name...) for each sample.

#### export

All captured data can be exported to CSV, PDF and HTML files for documentation.

#### virtual keyboard

To directly paste scanned barcodes into other applications, it is possible to activate the VirtualKeyboard to emulate keyboard input into the currently active application.





**QuickONE Export**

**Id:** 0087656643  
**Name:** Mayer, Markus  
**Birthday:** Freitag, 23. Februar 1979

Barcode	Type	TimeStamp
200000004131	DataMatrix	2014-11-14T11:48:49.2587
200000004055	DataMatrix	2014-11-14T11:48:54.2475
200000004049	DataMatrix	2014-11-14T11:48:55.0612
200000004058	DataMatrix	2014-11-14T11:48:57.6420
200000004051	DataMatrix	2014-11-14T11:49:02.4478
200000004139	DataMatrix	2014-11-14T11:49:03.1874
200000004137	DataMatrix	2014-11-14T11:49:05.7003
200000004138	DataMatrix	2014-11-14T11:49:06.5220
200000004134	DataMatrix	2014-11-14T11:49:09.0659
200000004141	DataMatrix	2014-11-14T11:49:09.0659
200000004133	DataMatrix	2014-11-14T11:49:09.7084

0087656643 - Editor

QuickONE Report

file:///D:/USER/MarcoLanzinger/Desktop/...

**QuickONE Export**

**Id:** 0087656643  
**Name:** Mayer, Markus  
**Birthday:** Freitag, 23. Februar 1979

Barcode;Type;TimeStamp  
200000004131;DataMatrix;2014-11-14T11:48:49.2587187+01:00;  
200000004055;DataMatrix;2014-11-14T11:48:54.2475956+01:00;  
200000004049;DataMatrix;2014-11-14T11:48:55.0612271+01:00;

**QuickONE Export**

generated / developed by www.labo-diagnostics.de - QuickONE

**ID:** 0087656643  
**Name:** Mayer, Markus  
**Birthday:** Freitag, 23. Februar 1979

Barcode	Type	Time Stamp
200000004131	DataMatrix	2014-11-14T11:48:49.2587187+01:00
200000004055	DataMatrix	2014-11-14T11:48:54.2475956+01:00
200000004049	DataMatrix	2014-11-14T11:48:55.0612271+01:00
200000004058	DataMatrix	2014-11-14T11:48:57.642031+01:00
200000004051	DataMatrix	2014-11-14T11:49:02.4478344+01:00
200000004139	DataMatrix	2014-11-14T11:49:03.1874019+01:00
200000004137	DataMatrix	2014-11-14T11:49:05.7003304+01:00
200000004138	DataMatrix	2014-11-14T11:49:06.5220367+01:00
200000004134	DataMatrix	2014-11-14T11:49:09.0659102+01:00
200000004141	DataMatrix	2014-11-14T11:49:09.0659102+01:00
200000004133	DataMatrix	2014-11-14T11:49:09.7084032+01:00

## Technical Specifications

### Data Collector

<b>Operating System</b>	Windows 8.1 Professional
<b>2D Barcode Reader</b>	Standard and inverted 2D DotMatrix Barcodes
<b>Export Formats</b>	PDF, CSV, TXT, HTML
<b>Data Exchange</b>	USB2.0
<b>Hardware</b>	10,1" Tablet
<b>Dimensions (total)</b>	260 x 182 x 10 mm
<b>Weight (total)</b>	0.65 kg

## Order Information

### QuickONE PAD

Item-Nr.	Description
<b>QO-ONEPAD</b>	mobile data collector for 2D barcodes with integrated optical unit

### QuickONE PAD Options

Item-Nr.	Description
<b>QO-PATIENT</b>	optional patient module to capture additional patient information and data for each sample.
<b>QO-DOCKING</b>	optional desktop docking station for charging and data exchange

## Our Distributor

